

Model 610 Unique Features and Benefits

The Model 610 is the first and only automated high level disinfection system that offers process control, documented quality assurance and ease of use. The following features and benefits identify the unique advantages of this model when compared against existing practices or other high level disinfection equipment.

| ADVANTAGES | FEATURES AND BENEFITS |
|---|--|
| High Level Disinfection Process Control and Documentation | <ul style="list-style-type: none"> • Step by step process documentation. • Quality Assurance Documentation - Expanded data collection complies with new JCAHO, AORN and Infection Control recommendations. • Improves staff accountability. |
| Sophisticated User Interface | <ul style="list-style-type: none"> • Intuitive control panel. • Message alerts inform users when to conduct routine cleaning, preventive maintenance, refill cleaning solution, order supplies, and possible equipment performance or conditional performance issues that may need to be corrected to assure that medical devices will be effectively cleaned and high level disinfected and are safe for next use • Department preference menu |
| Efficacy and Cleaning Effectiveness | <ul style="list-style-type: none"> • Multi-purpose trays are adaptable to accommodate different device sizes. • New design increases water movement through devices keeping them fully submerged, assuring effective cleaning. • External water supply nozzle. |
| Software Controlled Process Management and Safety Assurance | <ul style="list-style-type: none"> • Automatic lid lock prevents interruption of any cycle, resulting in a certifiable wash and HLD cycle each and every time. • Energy efficient high/low temperature management. • Optimizes user safety by preventing contact with hot water. • Automatic cleaning solution injection guarantees that proper cleaning protocols are met. • Automatic system clean. • Improves ease and convenience. |
| Multi-Cultural Interface | <ul style="list-style-type: none"> • Select between English, Spanish or French to maximize operational familiarity and confidence in multi-cultural settings. |
| Low Operational Costs and Maintenance | <ul style="list-style-type: none"> • Fully insulated cabinet and lid. • Simple PM Kit extends machine performance and life span. |

High Level Disinfection Systems for Healthcare

Cenorin, LLC • 6324 South 199th Place • Suite 107 • Kent, WA 98032
Phone: 253-395-2400 • Fax: 253-395-2650 • Toll Free: 1-800-426-1042
FDA Registered • CSA Certified • ISO 13485:2003